

# EMT—ENDORING® MARKING TABS

EMT

## FUNCTIONS

The EMT – Endoring Marking Tabs enable the clinician to easily organize the Endoring Foam inserts before a procedure. The three colors of the EMT help differentiate files of varying lengths or types when they are placed into a foam insert prior to autoclave processing. The user can track how many times a set of files have been used by easily removing one of the six numbered tabs located on the top of the EMT. The textured surface on the EMT can be marked with an indelible ink pen to indicate patient or file set-up information.



## TECHNICAL SUMMARY

### EMT:

The EMTs are molded from a thermoplastic resin that can withstand limited steam autoclaving. Each EMT is molded in three semi-transparent colors, including green, red and blue. The design allows for easy insertion into any Endoring Foam Insert prior to autoclave processing. A textured surface is provided for easy identification marking with a permanent marker such as a Sharpie® pen. The six numbered tabs are scored for easy removal by finger pressure or the use of an instrument such as a hemostat. Eight EMTs are attached to a convenient and unique dispensing wheel. Each EMT is easily detached from its wheel with a gentle twisting motion.

## PACKAGING

Each package contains three wheels of eight (24 total EMTs) in a sealed poly bag. Three dispensing wheels containing 8 EMTs are enclosed in each bag. All packages are labeled, CE marked and coded to HIBC convention.

## PRODUCT CODES

REF#	ITEM DESCRIPTION
EMT-s	Endoring Marking Tabs: Single package of 24 tabs (8 red, 8 green and 8 blue)
EMT	Case of Endoring Marking Tabs (24 bags/case)

ENDODONTIC ORGANIZING TOOL

FEATURES & BENEFITS

- Color Coded** Three colors are included in each bag to identify instrument set-ups. As an example, a blue EMT might mark 21mm files. Other colors included are green and red.
- Autoclavable** EMTs are molded from a resin that withstands limited steam autoclaving. Each EMT will last for one patient case or up to six autoclaving cycles. EMTs can also be processed with cold disinfectants.
- Easy Insertion** After placing files and other instruments into a foam insert, EMTs can easily be inserted to identify each set-up. The entire set-up can be bagged and placed into the steam autoclave for processing.
- Six Numbered Tabs** Six numbered tabs are located along the top edge of the EMT. Each tab can be snapped off with a hemostat or finger pressure to track each time the endodontic instruments are used.
- Convenient Marking Surface** The EMT's marking surface can be used to record patient names, dates, file sizes or any other important information required. Indelible ink pens, such as Sharpies, are ideal for this use.
- Unique Packaging** Eight EMTs are molded on a convenient dispensing wheel. There are three dispensing wheels per package (24 EMTs total).



800/752.2812  
(Inside the USA)

503/531.3904  
(Outside the USA)

503/531.3757 FAX

www.jordco.com

EMT



EMTs can be sterilized up to six times using Steam Autoclave. EMTs can also be placed in cold disinfectants.

Patent Pending  
Made in USA

