

PURE BOND® DISPENSER®



FUNCTIONS

Disposable, multi-purpose, light protected storage wells for bonding agents and composites. This efficient device can store light sensitive bonding agents, sealants and composites in four wells protected by a light locking lid. Bonding materials can be pre-loaded prior to a procedure and be easily delivered at chair side.

TECHNICAL SUMMARY

PURE BOND DISPENSER: Thermoformed black plastic. Approximately 2-1/4" by 2" with lid closed. 4 agent wells, 1 large utility well and 2 brush storage channels per unit.

PACKAGING

Supplied in a poly bag of 50 disposable units. All packaging is coded to HIBC convention.

PRODUCT CODES	REF#	ITEM DESCRIPTION
PBD-s		Pure Bond Dispenser: One package of 50 disposable dispensers
PBD		Case of Pure Bond Dispensers: 24 packages per case

BONDING AGENT & COMPOSITE DISPENSER



800/752.2812
(Inside the USA)

503/531.3904
(Outside the USA)

503/531.3757 FAX

www.jordco.com

COMPARISON GUIDE

FEATURES	PURE BOND	DANDY™	CENTREX®	GENERIC STORAGE WELLS
Number of agent wells	4	4	4	2 or more
Number of utility wells	1	0	0	0-?
Sealed wells to prevent spills	YES	NO	NO	NO
Light safe lid	YES Integral Light Lock lid	YES Separate loose fit	YES Separate loose fit	?
Hand held	YES	YES	YES Awkward	YES
Flat surface operation	YES One hand only	YES	YES Awkward	YES
Brush storage channels	YES: 2	YES	NO	NO
Disposable	YES	YES	NO Top must be cold soaked. Base=disposable	YES & NO
Disinfection required	NO	NO	YES: reusable top	YES: If not disposable
Non-reactive resin	YES Stable in all bonding systems	?	?	?
Compact & easy to handle	YES Lid attached	NO Loose lid	NO Loose Lid	VARIABLE

FEATURES & BENEFITS

- Four Agent Wells** Deep wells accommodate all agents needed for bonding procedures. Great for etchants, bonding agents, sealants and composites.
- Light Protecting Lid** Unique, one piece hinged design prevents any light contamination of bonding agents or composites. Light sensitive materials can be stored for several hours without activation.
- Leak Proof Lid** Manufactured of black thermoformed plastic. The locking lid prevents accidental spillage of bonding agents.
- Brush Storage Channels** Dispenser channels store small applicator brushes prior to a bonding procedure.
- Disposable, Single Use** Hygienic single use tray is disposable after each use. Reduces cross contamination and saves valuable chair and set-up time.
- Increased Efficiency** Assistants can pre-load required dispensers for the day. The light protecting lid will keep the sensitive bonding materials fresh. The original containers for the bonding agents do not have to be disinfected after each procedure.
- Compact and Easy Handling** Convenient size can be easily held by an assistant or placed on any surface. Easy open and close lid simplifies the procedure.
- Saves Bonding Agents** Unique well design pools remaining agents for optimum utilization.
- Durable Resin Construction** Compatible with most resins, primers, solvents and etchants. Prevents contamination of agents by dissolved resin precipitates.



Properly dispose of dispenser after use to help reduce cross contamination. Materials can be placed for extended storage with the lid closed. The light lock lid prevents any light from contaminating agents.

U.S. Patents 5,749,730
Made in USA

