

# PURE FLOSS DISPENSER®

SANITARY DUAL FLOSS DISPENSING SYSTEM



### FUNCTIONS

A unique floss dispensing system in which the floss is suspended between the storage and cutting towers. This suspension of floss allows the operator to use a no-touch technique to dispense floss and reduce cross contamination.

### TECHNICAL SUMMARY

- BASE:** Advanced non-skid elastomer (gray). Can be cleaned with most disinfectants.
- BODY:** Polymer coated heavy zinc casting (almond). Holds two spools of standard bulk floss. Can be cleaned with most disinfectants.
- CUTTER:** A removable, autoclavable, single piece, stainless steel cap with dual cutters.
- SIZE:** Approximately 4" x 6" footprint, 3" height.

### PACKAGING

Available as a single unit. Most offices will require one per operatory. All units packaged in a corrugated box and coded to HIBC convention.

### PRODUCT CODES

REF#	ITEM DESCRIPTION
PFD-s	Pure Floss Dispenser: 1 base/holder assembly and 1 removable cutter cap
PFD-cc	Pure Floss Dispenser Cutter Cap (replacement)
PFD	Case of Pure Floss Dispensers: 12 dispensers per case



800/752.2812  
(Inside the USA)

503/531.3904  
(Outside the USA)

503/531.3757 FAX

www.jordco.com

COMPARISON GUIDE

FEATURES	PURE FLOSS	DISPOSABLE BRAND	HAND-HELD PLASTIC BULK
No-touch dispensing	YES	NO	NO
Countertop stable	YES	NO	NO
Replaceable floss cutters for optimal performance	YES	N/A	NO
Removable and autoclavable floss cutters	YES	NO	NO
Dual floss dispensing	YES	NO	NO
Modern design eliminates contamination traps	YES	NO	NO
Injection molded, sealed non-skid base	YES	NO	NO
Smooth polymer coated finish for easy disinfection	YES	NO	NO
Fast floss dispensing	FASTEST	N/A	SLOWEST
Fast floss dispensing during procedures	FASTEST One hand, no-touch technique	SLOW Requires 2 hands	SLOW Requires 2 hands + contaminates housing
Reduces cross contamination	YES	YES	NO

FEATURES & BENEFITS

- Heavy Construction** The heavy weight of the body and skid-resistant base combine to make a very stable unit. The Pure Floss Dispenser will not move during normal dispensing, but is easily portable for cleaning purposes.
- Sanitary Design** Pure Floss Dispenser has a precise fit and finish which reduces areas that can trap contaminants. This elegant design reduces cross contamination while simplifying cleaning and disinfection.
- No-Touch Dispensing** The dental floss is suspended between the storage and cutting towers. The floss is picked up with one hand to free it from the cutter, the desired length of floss is pulled from the storage tower, and the floss is cut, leaving non-contaminated floss for the next user.
- Easy Bulk Floss Loading** The advanced elastomer base provides excellent stability and will not mark countertops.
- Durable Metal** Polymer-coated, die cast zinc base and precision stainless steel cutter are built for durability. This is the type of construction and finish used in only the best dental equipment available today.
- Industry's Choice** Recommended by: *The Dental Advisor* which rates the Pure Floss Dispenser "+++1/2" on a scale from 1-5. Also recommended by *Reality* which awarded the Pure Floss Dispenser "Reality's Choice" by receiving a rating of "above 4" on a scale from 1-5.



To disinfect the surface, use any acceptable surface disinfectant recognized by the EPA. However we do recommend that you use caution with iodophors as they may cause some slight discoloration. Apply the disinfectants with a soft cloth or a cotton sponge. You may remove the stainless steel cutter cap and autoclave it as needed.

U.S. Patent DES. 366,537  
Made in USA

